DIMA

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Important:

Please, read the Instructions for Use carefully.
 Not following the indications of the user instructions may cause malfunctioning of the system or cause damage to the patient.


This document pretends to explain the instructions for use of the Surelift Uplift. These instructions for use are not a document of reference to correct the vaginal prolapse. Surelift Uplift should only be used by doctors instructed in the treatment of this kind of pathology and more specifically in the use of this device and following the instructions for use.

These instructions for use explain the common steps to implant the system. The procedure could change when using the system in special techniques.

Intended use:

Surelift Uplift is a system for the treatment of the vaginal prolapse associated or not to other pathologies. The origin of this pathology is the weakness of the tissues that support these structures maintaining them in their position in the pelvic floor. The system reconstructs the support structures of the pelvic floor repositioning the organs in its natural position.

Description of the device:

<p>SURELIFT UPLIFT</p> <p>Ref: UPLIFT</p> <p>It is composed by a monofilament polypropylene mesh. When it stays implanted, the result is reinforcing the pelvic floor by inducing a fibrotic tissue through and around the mesh.</p>	
--	--

<p>SURELIFT UPLIFT KIT</p> <p>Ref: KITUPL-05</p> <p>It is composed by an anchoring handle, a monofilament polypropylene mesh and five anchors in a charger support. When it stays implanted, the result is reinforcing the pelvic floor by inducing a fibrotic tissue through and around the mesh. The function of the anchors is to fix the mesh to the promontory.</p>	
---	---

Description of the accessories used with the device:

Accessories - Not included - Sale separately:

<p>ANCHORING HANDLE</p> <p>Ref: APLIC05</p> <p>The Anchoring handle allows the easy placement of the Anchors. The Anchoring handle is provided sterile and is reusable.</p>	
---	--

<p>UPLIFT ANCHORS</p> <p>Ref: UPL-05A</p> <p>A charger support with 5 sterilized anchors. The anchors are single use.</p>	
---	--

<p>ANCHORING HANDLE KIT</p> <p>Ref: UPLAPLIC5A</p> <p>A sterilized Anchoring Handle and a charger support with 5 sterilized anchors. The Anchoring Handle and the anchors are single use.</p>	
---	--

Warnings and precautions:

- 1.- The device must only be used by doctors trained in the placement of the Surelift Uplift. It is mandatory to know the mechanism of action of the device.
- 2.- The device must be implanted only by doctors who are qualified to develop the required surgical technique.
- 3.- This product is presented sterile and must be sterile before its use. Before opening the product the package must be inspected to check it is not damaged and sterility is not compromised.
- 4.- This product is not reusable. Surelift Uplift cannot be re-sterilised. Reject and do not use any unpacked system.
- 5.- This product is an implant, so strict aseptic measures must be taken during the surgical procedure.
- 6.- Inform the patient that disuria, hematuria or any other problem should be communicated to the surgeon as soon as possible.
- 7.- Appropriate surgical procedures must be performed to implant Surelift Uplift as well as for the treatment of contaminated or infected wounds.
- 8.- Retro pubic post-operative bleeding may occur. Observe any signs or symptoms before discharging the patient.
- 9.- Cystoscopy should be performed to confirm bladder integrity after the needle-passage.
- 10.- It is recommended that patients do not do big efforts or effort-sports (for example biking, jogging, etc.) at least for a minimum of three or four weeks. It is also recommended avoiding sexual intercourse during the first month after surgery.
- 11.- The date of the last actualisation of this document is included at the end of the Instructions For USE. If more than 24 months have passed since the last revision date, please contact the manufacturer or distributor to receive the updated Instructions For Use.

Adverse reactions:

- 1.- Vessels, nerves, bladder, urethra or intestine punctions or erosions may occur during the pass of the needles. The repair of the damage may require surgical intervention.
- 2.- Temporary surrounding tissue irritation or foreign body response may happen. These reactions may produce extrusions, fistula formation and/or inflammation.
- 3.- As a foreign body, Surelift Uplift could boost a pre-existing infection.
- 4.- In the case the patient will suffer complications or reactions caused for any of the components, there will have to be explanted immediately.
- 5.- The device may produce pain, mucosal or other structures erosion and/or migration of the mesh that may result in the failure of the intervention, recurrent prolapse and / or incontinence.

Contraindications:

- 1.- Do not implant the device in patients who are undergoing anticoagulation treatments.
- 2.- Do not implant the device in patients who suffers urinary tract infections or obstructions.
- 3.- Do not implant the device in pregnant women or women that are considering future pregnancies.
- 4.- Do not place the device if the patient suffers renal insufficiency.
- 5.- Do not place the device when autoimmune sicknesses affect connective tissues.

Surgical Technique:

Before using Surelift Uplift the doctor must read and understand the instructions for use.

Anaesthetics and antibiotic therapy:

The surgical technique to implant **SURELIFT UPLIFT** can be performed under regional or general anaesthesia. It is recommended to follow the antibiotic protocol dictated by the hospital.

SURELIFT UPLIFT placement: (Laparoscopic approach)

- 1.- By laparoscopic technique, to achieve correct exposure of the pelvic organs, mobilize ascending and descending colon.
- 2.- If the uterus is preserved, fixation of the uterus to the abdominal wall through percutaneous suture with needle straight. Pass the suprapubic suture through the abdominal wall, pass through the uterine fundus and pass out again the suture through the abdominal wall. The suture is tied with security on a piece of gauze.
- 3.- Identify the sacral promontory and push the peritoneum on him, avoiding the nearby vessels.
- 4.- Continuous the incision of the lateral peritoneum rectum to expose the pelvic muscles of the fund of the pelvis. It must be preserve the peri-rectal fat to avoid damage. Repeat dissection on the opposite side.
- 5.- Set the extreme form of "Y" mesh UPLIFT bilaterally to the elevators of the anus. It is recommended the use of any method for the correct identification of the vaginal dome.
- 6.- Release the suture which moves the uterus to the abdominal wall and open a window in the broad ligament (right), avoiding the Uterine artery and Fallopian tube.
- 7.- Identify the boundary of the vaginal dome to guide the dissection of the bladder from the vagina. Dissecting until guess the ball of the probe.
- 8.- Pass the end loose mesh through the window of the broad ligament, and making sure that it is not coiled on itself in all its way, fix it to the area of former dissection stitched to the apex and on both sides.
- 9.- With a forceps pull from the central part of the mesh on the promontory applying the desired strength and set the mesh through two anchors or suture to the anterior vertebral ligament,.
- 10.- Cut the excess of mesh and close the peritoneum avoiding any mesh exposure.

Instructions to precharge the Anchoring Handle:

Step 1: Anchoring Handle initial position.

Step 2: Introduce the Anchoring Handle in the charger support, matching up the tube slots with the "+" axles of the anchor. Press the grip to fix the anchor in the Anchoring Handle.

Step 3: Get out the Anchoring Handle from the charger support. The anchor stay fixed to the Anchoring Handle.

**Materials used in manufacturing:**

Monofilament Polypropylene Mesh:
Implantable polypropylene monofilament mesh

anchors:
PEEK

Sterilization:

Surelift Uplift is sterilised by Ethylene Oxide.
This product is a single use device and it must not be re-sterilized or re-used.
Do not use if the package has been damaged.

Storage:

It is recommended to store the Surelift Uplift at a temperature lower than 40°C, far from humidity, heat and direct light. Do not use after the expiration date detailed in the package.

Single use product:

Device intended to be used once only for a single patient. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

Last Update of the Instructions for Use:

Date of the last review: October 2011

Manufacturer information:

Manufacturer:	Desarrollo e Investigación Médica Aragonesa S.L.
Manufacturer address:	Polígono Industrial Mediavega, Parcela 2.8 50300 Calatayud (Zaragoza) Spain Tel: +34-976-884065 Fax : +34-976-891872 e-mail: dimasl@dimasl.com

Distributor information:

Distributor:	Neomedic International S.L.
Distributor address:	C/ Maestrat, nº 41-43 1º 08225 Terrassa (Barcelona) Spain Tel: +34-93-7804505 Fax: +34-93-7806105 e-mail: info@neomedicinternational.com