


AGERS® system


The Minimally Invasive **Vaginal Agenesis** System

THE MINIMALLY INVASIVE TREATMENT FOR THE ROKITANSKY SYNDROME



- Natural Neovagina formation
- Aprox.1 cms longer canal per day
- Comfortable adjustment by the patient

Agers® Vaginal Agenesis System



- **AGERS is a pulling mechanical system to invaginate the membrane creating a new vagina for those patients with a very short or undeveloped vagina (M.R.K.H. Syndrome*).**
- **It's placed by laparoscopy.**
- **No tissue grafts or complicated surgical interventions required.**
- **Traction is increased from the abdominal site, and manually adjusted by the patient every day.**
- **In 8 days may create a neovaginal canal up to approximately 7-8 cms.**

(*) **MRKH:** Mayer-Rokitansky-Kuster-Hauser syndrome

CODE	DESCRIPTION	INCLUDES
AGERS	Agenesis Remeex System	1 varitensor + 1 base support + 1 traction thread + vaginal balls + vaginal needle system MAGEO1

* Consult us for vaginal dilators.

info@neomedicinternational.com

© 2013 Neomedic International. All rights reserved.

www.neomedic.com

neomedic
international

The Continece Company